### New data analysis and development of new GRNsight functionality will aid researchers in discovering more about gene regulatory networks

Blair Hamilton, Department of Mathematics Emma Tyrnauer, Department of Biology Katie Wright, Department of Chemistry and Biochemistry Zach Van Ysseldyke, Department of AIMS

December 12, 2017 BIOL/CMSI 367-01: Biological Databases, Loyola Marymount University

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.


# Organisms have mechanisms that allow them to withstand environmental stresses

- Organism -- Shivering
- Organ Systems -- Endocrine, Circulatory
- Organ -- Digestion
- Tissue -- Muscle movement
- Cellular -- Cold shock proteins



#### Yeast is an ideal model organism for learning about Eukaryotic cells

- Significant growth in a matter of days
- Easy to manipulate in lab
- Can resist and respond to many different environmental stressors
- Controls cell division in a way similar to ours




### The right click functionality will be very useful for gene research

- Current Functionality: Visualize gene regulatory networks
- Goal: ability to easily open gene information page
- Solution: Adding right click functionality
  - Readily available gene information
  - No need for extensive database searching

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

#### Individual milestones contribute to group effort



- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.


Percentage of genes demonstrating a cold shock response decreased with stricter p-values

- Standard p < 0.05 revealed ~41% of genes demonstrating significant response to cold shock
- B & H corrected p < 0.05 revealed second largest percentage
- Bonferroni-corrected p < 0.05 was the most strict


ANOVA	WT
p < 0.05	2528 (40.85%)
p < 0.01	1652 (26.70%)
p < 0.001	919 (14.85%)
p < 0.0001	496 (8.01%)
Benjamini & Hochberg- corrected p < 0.05	1822 (29.44%)
Bonferroni-corrected p < 0.05	248 (4.01%)

### Genes were clustered into groups with similar regulation mechanisms

- Clusters ordered by significance
- Cluster 45 and 48 genes demonstrate upregulation
- Cluster 9 and 0 demonstrate downregulation



#### Profile 45 demonstrates clear upregulation after cold shock



### A table of GO terms sheds light on the function of the regulatory genes


Category Name	#Genes Category	#Genes Assigned	#Genes Expected	#Genes Enriched	p-value	Corrected p-value	Fold
cellular component organization (GO:0016043)	540	161	162.7	-1.7	0.6	1	1
nuclear lumen (GO:0031981)	327	183	98.5	84.5	2.40E-27	<0.001	1.9
<b>RNA binding</b> (GO:0003723)	215	116	64.8	51.2	6.30E-15	<0.001	1.8

### Transcription factors for the network were selected based on significance

TF	NDT80	YAP1	PDR3	UME6	MIG2	ZAP1	HAP4	CIN5
P-value	9.74E-10	1.31E- 09	2.82E- 09	6.12E- 09	1.96E- 08	2.99E-0 3	3.02E- 02	0.806


TF	YOX1	SFP1	YHP1	SUT1	STB5	ACE1	MSN2	GLN3
P-value	0	0	0	7E-15	7E-15	1.9E-14	2.72E-13	9.15E-11

### Some genes are interconnected more than others



Arrows indicate interaction present

### A weighted network details the nature of the interaction



- Pink lines indicate induction
- Blue lines
   indicate
   repression
- Thickness of lines indicate magnitude of interaction

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

### Right click functionality allows for easy access to additional gene information

```
.on("dblclick", nodeTextDblclick)
.on("contextmenu", function (gene) {
   var tempLink = $("<a></a>")
        .attr({
            href: "/gene/info.html?" + $.param({symbol: gene.name}),
            target: "_blank"
        });
   $("body").append(tempLink);
    tempLink.get(0).click();
    tempLink.remove();
   d3.event.preventDefault();
```



### Inter team collaboration aided in development of new functionality

Web app and service for visualizing models of gene regulatory networks

S000005111 etc. YNL167C <name xmIns="http://uniprot.org/uniprot">SKO1\_YEAST</name> MA0382.1

#### **SKO1Genus Species**

GRNSight

General Information
DNA Sequence
Protein Information
Gene Map
Regulation
Interaction
Gene Ontology
Frequency Matrix and Sequence Logo

### Inter team collaboration aided in development of new functionality

var sgdHrefTemplate = "https://www.yeastgenome.org/locus/"; var sgdId = gene.sgd.sgdID; \$(".sgd-link").text(sgdId).attr({ href: sgdHrefTemplate + sgdId });

// DNA Sequence Tab
var ensemblDNA = gene.ensembl.dnaSequence;
\$(".dnaSequence").text(ensemblDNA).attr({ href: ensemblHrefTemplate + ensemblDNA });

#### Frequency Matrix and Sequence Logo help visualize data from JASPAR



 A:0
 0
 0
 4921

 C:0
 0
 100841518

 G:1001000
 163728

 T:0
 0
 0
 0
 32

### Code for Frequency Matrix and Sequence Logo

var frequencyMatrix = gene.jaspar.frequencyMatrix; var a = ""; for (var i = 0; i < frequencyMatrix.A.length; i++) { a += "" + frequencyMatrix.A[i] + ""; }

var sequenceLogo = gene.jaspar.sequenceLogo;
\$(".sequenceLogo").attr({ src : sequenceLogo });

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

# This team is confident that our contributions will aid future research

- Studying wild type yeast can show importance of certain genes in cold shock response when compared to mutated cells.
- Dynamic integration enhances user interface by providing the user with in depth information on the desired gene.
- Newly analyzed data and newly coded functionality will both lead to new discoveries about gene regulatory networks.

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

### There are many avenues available for further research

- Microarray data of cold shock response in other eukaryotes
- Microarray data of other stress responses in yeast
- Prediction of related genes once right click function has been executed
- A function on the new gene page to alert the user that the page is loading or that data is being pulled

#### Summary

- Coding new GRNsight functionality and performing yeast data analysis can lead to new insights in gene regulatory research.
- Achieving Individual milestones contributes to group effort.
- Cold shock regulation differs between gene clusters.
- Right click function links to a new gene page which helps biologists find information quickly.
- This team's work will aid future research
- There are many avenues available for further research.

#### Acknowledgements

Loyola Marymount University, Seaver College of Science and Engineering

The Departments of Biology, Mathematics, Computer Science, Chemistry and Biochemistry, and AIMS

**Professor Dionisio and Professor Dahlquist** 

**The Eddies** 

Students of BIOL/CMSI 367-01: Biological Databases



#### References

Get url parameter jquery Or How to Get Query String Values In js. (n.d.). Retrieved December 04, 2017, from https://stackoverflow.com/questions/19491336/get-url-parameter-jquery-or-how-to-get-querystring-values-in-js

How to retrieve GET parameters from javascript? (n.d.). Retrieved December 04, 2017, from

https://stackoverflow.com/questions/5448545/how-to-retrieve-get-parameters-from-javascript JavaScript Window Location. (n.d.). Retrieved December 04, 2017, from

https://www.w3schools.com/js/js\_window\_location.asp

Murata, Y., Homma, T., Kitagawa, E., Momose, Y., Sato, M. S., Odani, M., Shimizu, H., Hasegawa-Mizusawa, M., Matsumoto, R., Mizukami, S., Fujita, K., Parveen, M., Komatsu, Y., Iwahashi, H. (2006). Genome-wide expression analysis of yeast response during exposure to 4 C. Extremophiles, 10(2), 117-128.
Window.location. (n.d.). Retrieved December 04, 2017, from

https://developer.mozilla.org/en-US/docs/Web/API/Window/location